


**POC SOLUTIONS**  
**OPTIMISE PATIENT CARE**  
With semi quantitative results


**HCG - FOB - PSA - CRP**

**For professional use**

**NEW PLATFORM**



**INCREASED DETECTION RANGE**



Negative  
Low positive  
Positive  
High positive




## SEMI QUANTITATIVE TESTS

NG-semiQ tests	Specimen	Format	Storage	Shelf life	Cat. Ref.
NG-semiQ hCG*	Urine/Serum	Device	4-30°C	24 months	NGS-HCG-C23
NG-semiQ FOB*	Feces	Device	4-30°C	24 months	NGS-FOB-F23
NG-semiQ PSA*	WB/S/P	Device	4-30°C	24 months	NGS-PSA-W23
NG-semiQ CRP*	WB/S/P	Device	4-30°C	24 months	NGS-CRP-W23

\* under development

## GENERAL TEST PROCEDURE



## INTERPRETATION

