

OTC Pregnancy Test

WITH CONCEPTION INDICATOR

It tells me more !



PRODUCTS

Product	Specimen	Format	Storage	Shelf-life	Cat. Ref.
NG-Test hCG OTC*	Urine	Midstream	4-30°C	24 months	NGS-hCG-U31-001

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity :


The sensitivity of the NG Biotech OTC pregnancy test (Urine) is 10 mIU/mL.

The specificity of the NG Biotech OTC pregnancy test is determined by measuring the cross reactions with the luteinising hormone (LH) and the follicle stimulating hormone (FSH). With respective concentrations of 200 mIU/mL of LH and 200 mIU/mL of FSH, the results were negative.


Standardisation :

The NG Biotech OTC pregnancy test (Urine) is standardized against NIBSC hCG reference standard (Ref: 07/364). Additional assessments with human urine samples with known hCG concentrations were performed showing a 99.99% correlation with usual laboratory standard techniques as well as inter/intra batch reproducibility.

GENERAL TEST PROCEDURE



INTERPRETATION



*Coming soon